

A validated solution for shipping primary human proximal tubule cell monolayers between laboratories



Kidney

Lucy Gentles, Claire Devlin, Cindy Carr, Colin Brown

Newcells Biotech, Newcastle Helix, Newcastle Upon Tyne


Poster number: P91

Introduction

- aProximate™ is a highly sophisticated, near-physiological proximal tubule cell (PTC) model that is a proven market-leader in the study of transporters and drug-drug interactions.
- One issue with primary cell cultures is the short time window they remain differentiated and functional, which has hindered the sharing of primary cell models between laboratories.
- To overcome this and ship cells globally in an assay-ready format, a safe and reliable short-term preservation method was developed and internally validated, pausing cells at hypothermic temperatures for up to 5 days.

Methods

- Human PTCs were isolated from the cortex tissue of 17 nephrectomies by Percoll gradient separation and grown on Transwell® filter supports to functional maturity (Figure 1.)
- Membrane barrier integrity was quantified by transepithelial electrical resistance (TEER), with values of more than $60~\Omega.cm^2$ used as quality control to indicate a fully established, polarized monolayer. Cultures typically took 5-6 days to reach confluence.
- Cells were placed into metabolic stasis, using a formulated biogel and maintained at 4°C for 3 or 5 days to simulate shipping conditions.
- Subsequently, cells were revived to 37°C and left to recover until the minimum TEER QC value was reached.



<u>Figure 1.</u> *Illustration of methodology* Cells were isolated by enzymatic digestion and Percoll gradient preparation, cultured to confluence, then coated in a formulated biogel and placed into hypothermic preservation for up to 5 days. Upon revival and TEER recovery, monolayer viability and functionality were assessed.




Figure 2. Membrane barrier integrity assessment. (A) TEER recovery was cold-time dependent, with plates shipped for 3 (n = 7 kidneys) and 5 days (n = 12 kidneys) taking 24hr and up to 72hr respectively to reach QC threshold. (B) Lucifer Yellow (LY) paracellular leak in monolayers shipped for 5 days was 0.32 % (n = 4 kidneys) and remained well below the upper accepted limit of 2%.





Figure 3. Morphology & Protein Expression. (Phase Contrast) Following 5 days of hypothermic preservation and a 48-72hr recovery phase, monolayers remained differentiated and exhibited a single-layer cobblestoned epithelium, morphologically akin to pre-shipment control. (ZO1) Extensive positive immunofluorescent staining of the tight junction protein, ZO1 (red), was demonstrative of an uncompromised cell membrane barrier. (Annexin V) Expression of the apoptotic marker, Annexin V was not affected by shipment. (SLC22A6 / OAT1) & (SLC22A8 / OAT3) Expression of Organic Anion Transporters 1 and 3 were not affected by shipment. All nuclei counterstained with Hoescht (blue) and images taken on the ImageXpess Pico (Molecular Devices).






Figure 4. (A) Cell health assessment. Monolayers were viable and metabolically active, indicated by way of ATP assay (>100 % of control, n = 9 kidneys). LDH cytoxicity revealed cell loss of approximately 26% upon initial monolayer revival but PTCs grew back to full confluency during the 1-3 day recovery phase. (B) Transporter Expression. At the mRNA level, expression of the key drug transporters OCT2, MDR1, MATE1 and MATE2K were stable or upregulated relative to pre-shipment control (93%, 117% 123% and 128% respectively, n = 5 kidneys).




Figure 5. Bidirectional transcellular flux and uptake of a radiolabelled compound (A) Functionally, monolayers conserved the mechanics to facilitate the OCT2-mediated basolateral uptake and MATE-mediated efflux of 14 C creatinine. The efflux ratio for creatinine in shipped cells was 1.8 ± 0.16 fold and was not significantly different from time-matched control monolayers (n = 10 kidneys, paired t test, P > 0.01). (B) Likewise, the basolateral: apical creatinine uptake ratio was 3.5 ± 0.43 and was not significant from time-matched control monolayers (n = 17 kidneys, paired t test, P > 0.05).

Conclusion

Taken together, these data suggest that coating PTCs at hypothermic temperatures offers the possibility of shipping a viable and functional primary cell model across the world with improved confidence, overcoming a major limitation in the widespread use of primary cell models.