
Comparison of multiple batches of human iPSC-derived retinal organoids produced at large scale

Purpose

Due to well-known limitations of in vivo and existing in vitro retinal models, a 3D in vitro model of the human retina which is reproducible and able to accurately predict in vivo outcomes is highly desirable.

Our aim consistency of human iPSC-derived retinal organoids (RO) produced at large scale by quantifying the and gene expression levels of key retinal cell markers across differentiation in multiple batches.



Results




Figure 1. Photoreceptor (PR; Genes: RCVRN, RBP3, IMPG1, CRX; Protein: RCVRN), bipolar (BC; Genes: GRIK1, CADPS; Protein: PKCα), Müller glia (MG; Genes: RLBP1, CRYM; Protein: CRALBP), retinal ganglion cell (RGC; Genes: MATH5, BRN3; Protein: SNCG), horizontal and amacrine (HC/AC; Genes: PROX1, TFA2A; Protein: PROX1, AP2α), Cone PR (Genes: OPN1SW, OPN1MW, OPN1LW, ARR3, RXRG; Protein: OPN1MW/LW) and Rod PR (Gene: RHO, NRL; Protein: RHO) cell markers are expressed at different times throughout RO development which resembles in vivo development. (A) Gene expression and (B) protein expression of RGC markers peak between day 30-90, when PR progenitor cells start differentiating. Expression of cone and rod PRs initiate from day 120 and achieve highest expression and maturation levels after 210 days in culture.

Carolina Gándara¹, Will Atkin¹, Hannah Steward¹, Madalaine Kay¹, George Buchanan¹, Lyle Armstrong², Valeria Chichagova¹ 1. Newcells Biotech, Newcastle upon Tyne, UK. 2. Newcastle University, Newcastle upon Tyne, UK

Conclusions

We analysed the gene and protein expression profile of key retinal cell markers across differentiation in four batches of human iPSC-derived ROs. We observed that PRs, HCs and ACs, at later stages of development, and RGC at early stages, appear at consistent levels across multiple batches. This data set provides crucial information for pre-clinical studies in ROs with application in drug discovery, disease modelling and gene therapy.

