

Development of a predictive in vitro human proximal tubule platform for de-risking nephrotoxicity Piyush Bajaj², Git Chung¹, Keith Pye¹, Mike Nicholds¹, Yvonne Dragan³, Matthew Wagoner², Colin Brown¹ ¹Newcells Biotech Ltd Newcastle upon Tyne UK, ²Global Investigative Toxicology, Drug Safety Research & Development,³ Discovery Toxicology, Drug Safety Research & Development Takeda Pharmaceuticals, Cambridge, MA USA **Toolkit of 36 Compounds Tested**

Introduction

- Nephrotoxicity is a major reason for drugs failing during clinical development.
- Currently there is no in vitro platform that enables cross-species comparisons of drug transport or nephrotoxicity.
- Our innovative solution is to develop highly differentiated assay platforms using primary renal proximal tubule cells(PTCs) derived from key animal species to measure both drug transport drug induced kidney injury a range of biomarkers across species
- we showcase a data from our highly • Here differentiated Human Primate proximal tubule model

Methods

- PTCs were isolated from fresh Human kidneys and cultured onto Transwell inserts as outlined below
- Test toolkit was made up of 36 compounds with known nephrotoxicity liability and clinical data were screened using the model Toolkit was made up of 19 Primary PT-toxic and 17 Secondary PTtoxic (3) or non-nephrotoxic (14) compounds
- Monolayer were exposed for 72hours to a range of concentrations (0-300µM) of test compound
- Toxicity was measured used 6 parameters of cell Resistance (TEER), Transepithelial health: Intracellular ATP concentration, LDH release , KIM-1 release, Clusterin release NGAL release









Diverse Chemical Structure

Physiochemical properties of Toolkit of Compounds Previously published dataset Takeda dataset **Range of Mechanistic Exposure Routes** Non - Nephrotoxic Non - Nephrotoxic Nephrotoxic 500 500₁ 6:3 (2 300 45 200 104 100 2:6 -12



Varma et al Pharm Res. 32:3785-3802 (2015)

Toolkit was chosen to encompass a diverse range of chemical structures and be substrates for a wide range of proximal tubule transporters

lomerulus	Cystatin C		
	KIM 1		
	Clusterin		
imal Tubule	NGAL		
	NAG		
	Osteopontin		



August 2018- FDA approved a safety composite biomarker panel to aid in the detection of kidney tubular injury in phase 1 trials in healthy human volunteers. Biomarkers are detected in the urine

Biomarkers were 3 from 5 qualified biomarkers of proximal tubule damage outlined by FDA

Biomarker Response to Nephrotoxic Challenge is Highly Predictive of In Vivo Outcome

Human proximal tubule cell monolayers retain a remarkable degree of differentiation and express a range of functional transporters and clinically relevant biomarkers of nephrotoxicity that are sensitive to nephrotoxin challenge over time. Human PTC monolayers show excellent potential as an in vitro predictive screening platform

cLogP + tPSA 89.5 42.9 33 Compounds were selected to minimize influence of physiochemical properties in

cLogP

57.9

73.7

Sensitivity (%) Specificity (%)

64.3

64.3

weighting the outcome of assay

cLogP < 0.8

tPSA > 104

Conclusions



nly AC ₅₀ value							
-	ATP	KIM-1	NGAL	Clusterin			
)	< 50	< 30	< 22	< 95			
6	0.61	0.61	0.64	0.66			
3	47.4	47.4	47.4	63.2			
2	88.2	88.2	88.2	88.2			

In vitro safety margin (AC₅₀/C_{max})

ł	ATP	Clusterin	KIM-1	NGAL
0	< 25	< 31	< 30	< 30
2	0.64	0.75	0.74	0.78
5	42.1	57.9	63.2	63.2
2	88.2	88.2	88.2	88.2

Injury-specific biomarkers (KIM-1, NGAL, Clusterin) showed better predictivity

Incorporation of exposure (total C_{max}) improved predictive performance of most