CHARACTERIZATION OF CISPLATIN TOXICITY IN aPROXIMATE™ HUMAN PROXIMAL TUBULE CELL MONOLAYERS




Git Chung, Lyle Armstrong, Mike Nicholds, Colin Brown* **Newcells Biotech Ltd UK and *Newcastle University UK**




Introduction

- Around 50 % of preclinical toxicity screens fail to predict subsequent toxicity in vivo.
- This leads to significant attrition of drug molecules during drug development.
- Understanding nephrotoxicity has been hampered by the lack of good renal model.
- Here we demonstrate the utility of the



Methods

- PTCs were isolated from human kidneys as described in Figure 1, and cultured onto Transwell inserts.
- Confluent monolayers were treated with nephrotoxins.
- Transepithelial electrical resistance (TEERs) were measured, in addition to cell viability using an Real time



aProximate human proximal tubule cell (PTC) monolayers as an *in vitro* tool to investigate nephrotoxicity using clinically relevant biomarkers NGAL, KIM-1 and clusterin.

via density gradients

Tubular cell retrieved and cultured on Transwell inserts

Figure 1: PTC isolation and culture.

Glo^{™™} assay lacksquare

Culture media was collected from basolateral both apical and compartments and analysed for biomarkers using an ELISA approach.

Basolateral chamb


Figure 2: PTCs cultured on Transwell insert to recapitulate proximal tubule epithelial formation.

Biomarkers are released in to the apical media in a concentration dependent manner

Apical medium

Cell lysate

Basolateral medium



cisplatin.

Apical medium Cell lysate 50 μM 10 μM 0 μΜ 1 μM 100 μM Cisplatin Cisplatin Cisplatin Cisplatin Figure 4: NGAL levels by aProximate monolayers after 24 hour treatment with

cisplatin.

μ<u></u> 30-100 μM 10 μM 50 μM 1 μM 0 μΜ Cisplatin Cisplatin Cisplatin Cisplatin Cisplatin Figure 5: Clusterin levels by aProximate monolayers after 24 hour treatment with cisplatin.

- Levels of all 3 clinically relevant biomarkers KIM-1, NGAL, and clusterin increased within response to challenge with a range cisplatin concentrations.
- Importantly, as in vivo, biomarker release was predominately across the apical membrane than across the basolateral membrane (Figure 3, 4 & 5).

Time-dependent increase in biomarker release in response to cisplatin challenge

ng/ml

Control




Figure 6: KIM-1 levels by aProximate monolayers at different periods of exposure to 10 μM cisplatin.






Figure 7: NGAL levels by aProximate monolayers at different periods of exposure to 10 μ M cisplatin.



monolayers at different periods of exposure to 10 µM cisplatin.

- Levels of KIM-1, NGAL and clusterin in control treated monolayers increased over a and period of 5 days (Figure 6, 7 & 8).
- Increase in KIM-1, NGAL and clusterin with time demonstrate the cumulative effect of exposure to nephrotoxin challenge.

Response of Apoptosis Signalling Pathways in Human PTCs to Cisplatin Challenge



Gene symbol	Description	Fold regulation	p- value
BAG3	BCL2-associated athanogene 3	2.1	0.0007
BAK1	BCL2-antagonist/killer 1	4.1	0.02
BAX	BCL2-associated X protein	3.9	0.0002
BCL2A1	BCL2-related protein A1	11.1	0.02
BCL2L1	BCL2-like 1	3.2	0.001
BID	BH3 interacting domain death agonist	2	0.02
ВІК	BCL2-interacting killer (apoptosis-inducing)	5.6	0.03
BIRC3	Baculoviral IAP repeat containing 3	2.5	0.0001
	Caspase 14, apoptosis-related cysteine		
CASP14	peptidase	21.8	0.0003
	Caspase 7, apoptosis-related cysteine		
CASP7	peptidase	2.1	0.0009
CD70	CD70 molecule	10.7	0.02
CFLAR	CASP8 and FADD-like apoptosis regulator	2.2	0.0001
DIABLO	Diablo, IAP-binding mitochondrial protein	2.4	0.001
FAS	Fas (TNF receptor superfamily, member 6)	3.9	0.00003
	Growth arrest and DNA-damage-inducible,		
GADD45A	alpha	3.2	0.0003
	Lymphotoxin beta receptor (TNFR superfamily,		

Figure 9: mRNA levels of apoptotic pathway proteins in aProximate monolayers at 72 hours exposure to 10 µM cisplatin. Cisplatin challenge resulted in significant increases in genes associated with apoptosis at the qPCR level. These included the death receptors; Fas (3.9 fold) CD70 (11 fold) and TNFR (3.8 fold), intrinsic pathway mediators BAX (4.1 fold) together with AIF (1.9 fold) and DIABLO (2.35 fold) and ER stress pathway markers resulting in significant increases in downstream caspase production leading to apoptosis. Importantly cisplatin also stimulated mRNA levels of the anti-apoptotic proteins; caspase 14 (21.8 fold) and BCL2A1 (11.1 fold).

Conclusion

In summary, aProximate[™] human proximal tubule cell monolayers retain a remarkable degree of differentiation, express clinically relevant biomarkers of nephrotoxicity and signalling pathways that translate to an appropriate response to cisplatin challenge aProximateTM human proximal tubule cell monolayers show excellent potential as an in vitro predictive human toxicology screening platform during the drug development process