




RETINAL ORGANIDS



RETINA

An *in vitro*, light responsive retinal model for accurate predictions of *in vivo* outcomes

The Retinal Organoids model from Newcells recapitulates the complex structure of the retina

Product Specification

Format	10 organoids per 5ml microfuge tube 150 ml of optimized cell culture medium per 100 organoids 2 x 96 well plates per 100 organoids 2 x Pasteur pipettes
Cell Types	Organoids Rod and Cone photoreceptors Retinal ganglion cells (RGCs) Bipolar cells Horizontal cells Amacrine cells Müller glial cells RPE (expected availability 2022) Retinal pigment epithelial cells
Species	Human Rat (In development) Non-human primate (In development)
Available analytical readouts	Immunofluorescence analyses Gene expression by RT-qPCR Transcriptomic analysis by single-cell RNA sequencing Cytotoxicity assays Cytokine release Flow cytometry Electron microscopy Custom assays
Origin	Healthy donor Patient samples
Assay Window	≥30 days

- ✓ The organoids are ~1.3 mm in diameter and contain ~40,000 cells.
- ✓ Tested for different applications including gene therapy, toxicology and retinal disease modelling.
- ✓ They form primitive photoreceptor outer segments leading to responsiveness to light.
- ✓ They respond to known toxins similar to that seen *in vivo*.

Gene expression in retinal organoids through different stages

Cell type	Gene	Timepoint of appearance	Timepoint of peak expression
Retinal ganglion cells	<i>MATH5 (ATOH7)</i>	d30 - d180	d60
	<i>BRN3 (POU4F2)</i>	d30 - d210	d60
Horizontal and amacrine cells	<i>TFAP2A</i>	d30 - d210	d150
	<i>PROX1</i>	d30 - d210	d150
Bipolar cells	<i>GRIK1</i>	d30 - d210	d150
	<i>CADPS</i>	d30 - d210	d150
Photoreceptors	<i>RCVRN</i>	d60 - d210	d210
	<i>RBP3</i>	d60 - d210	d210
	<i>IMPG1</i>	d120 - d210	d210
	<i>CRX</i>	d60 - d210	d210
Cone photoreceptors	<i>OPN1SW</i>	d120 - d210	d210
	<i>OPN1MW</i>	d150 - d210	d210
	<i>OPN1LW</i>	d120 - d210	d210
	<i>ARR3</i>	d60 - d210	d180
Rod photoreceptors	<i>RXRG</i>	d60 - d210	d150
	<i>RHO</i>	d120 - d210	d210
RPE	<i>NRL</i>	d90 - d210	d180
	<i>RPE65</i>	d60 - d210	d210
Müller glia	<i>RLBP1</i>	d90 - d210	d210
	<i>CRYM</i>	d30 - d210	d210

Retinal organoid differentiation follows the developmental timeline of embryonic development of the retina with various cell types arising at different times in a sequential manner. Long-term cell survival is hampered by the limitations of the *in vitro* culture conditions and the limitations of the model itself. (e.g. absence of visual cortex).

Expected protein expression in retinal organoids at d150 and d180


Cell type	Cell marker	Protein expression at d150	Protein expression at d180	Protein localization at d150-d180
Photoreceptors	RCVRN	✓	✓	ONL
Retinal ganglion cells	SNCG and HuC/D	✓	✓	INL/GCL
Cone photoreceptors	OPN1MW/LW	✓ *	✓	ONL
Rod photoreceptors	RHO	✓ *	✓	ONL
Bipolar cells	PKC-α	**	✓	INL
Amacrine cells	AP-2α	✓	✓	INL
Horizontal cells	PROX1	✓	✓	INL
Müller glia	CRALBP	✓	✓	All layers

* Small number of developing rods and cones

** Expressed at transcriptional level




RETINAL ORGANOID




RETINA


Characterization of iPSC-derived Retinal Organoids




Localization and distribution of photoreceptors (RCVRN, green) and retinal ganglion cells (SNCG, red) in retinal organoids at d150.




Representative IF image of retinal organoid at d210. (A) Localisation and distribution of cone photoreceptors (OPN1MW/LW, green) and amacrine cells (AP-2α, red); (B) detailed view of image A, showing cone photoreceptors with primitive outer segments.



Localization and distribution of Müller glia cells (CRALBP, red) in retinal organoids at d180.



Localization and distribution of amacrine cells (AP-2α, red) in retinal organoids at d150.



Localization and distribution of horizontal cells (PROX1, green) in retinal organoids at d150

Resources

- Application of organoid technology for retinal disease modelling and drug discovery; Chichagova, Drug Target Review, June 2020
- [Application of organoid technology for retinal disease modelling and drug discovery \(drugtargetreview.com\)](https://www.drugtargetreview.com)
- Human iPSCs generate light responsive retinal organoids with variable and nutrient dependent efficiency; Hallam et al, Stem Cells, 2018, 36(10), 1535-1551
- [Human-Induced Pluripotent Stem Cells Generate Light Responsive Retinal Organoids with Variable and Nutrient-Dependent Efficiency - PubMed \(nih.gov\)](https://pubmed.ncbi.nlm.nih.gov/)
- Room temperature shipment does not affect the biological activity of iPSC derived retinal organoids; Georgiou et al, PLOS One, 15(6), e0233860
- [Room temperature shipment does not affect the biological activity of pluripotent stem cell-derived retinal organoids \(plos.org\)](https://plos.org/)
- Enhancing immune function of hiPSC derived retinal organoids by incorporating microglial cells; Chichagova et al, Investigative Ophthalmology and VISUAL Science, 2020, 61(7)
- [Enhancing immune function of hiPSC-derived retinal organoids by incorporating microglial cells | IOVS | ARVO Journals](https://iovs.arvojournals.org/)
- Human iPSC differentiation to retinal organoids in response to IGF1 and BMP4 activation is line and method dependent; Chichagova et al, Stem Cells, 2020, 38(2), 195-201
- [Human iPSC differentiation to retinal organoids in response to IGF1 and BMP4 activation is line- and method-dependent - PubMed \(nih.gov\)](https://pubmed.ncbi.nlm.nih.gov/)

The need for pre-clinical Retinal Organoids	Results with Newcells Retinal Organoids
Recapitulation of the complex architecture with the relevant cell types	A physiologically relevant, functional, light responsive model for mechanistic insights
Unlimited material for use in safety and efficacy studies	Simple pre-clinical studies for safety and efficacy
Lack of suitable models for the human retina for disease modelling for accurate pre-clinical data	Predictive disease modelling platform
Ethical responsibility related to the 3R principles (Replace, Reduce and Refine)	Reduction of use of animal models in line with NC3Rs and NA3RsC