

Utilization of high-content imaging for the study of lung fibrosis

Colin D.A Brown, Fiona Leslie, Oliver Birch, Chloe Whiting, Megan Webster


Newcells Biotech Ltd, The Biosphere, Draymans Way, Newcastle Helix, Newcastle upon Tyne. UK NE4 5BX



FMT is a pathologic mechanism of fibrosis


- Inhalation of toxic microparticles can damage lower airway epithelial cells. Sustained microinjury of the epithelium results in aberrant epithelial-fibroblast communication via the release of fibrotic stimuli, including TGF- β 1.
- TGF- β 1 drives the phenotypic transition of fibroblasts to myofibroblasts (FMT). Activated myofibroblasts are responsible for the excessive deposition of extracellular matrix (ECM) proteins, a characteristic feature of fibrosis. This leads to irreversible damage to the lung architecture and function. Therefore, new anti-fibrotic therapies are needed to suppress disease progression.

Modelling FMT processes *in-vitro* using high-content imaging enables compound screening to assist anti-fibrotic drug-discovery.




Optimization of culture conditions for assay robustness

Culture plate coatings were compared to assess fibroblast adherence and reduce assay variation




Significant cell detachment following mechanical plate washing was observed for coating 1 whilst coating 2 improved assay robustness.

Different media conditions were compared to achieve a substantial and consistent assay window



Media 3 promotes FMT, as indicated by the increased expression and quantification of α -SMA and collagen I following stimulation with TGF- β 1.

Extracellular Collagen I




The inclusion of a macromolecular crowder (MMC) in the assay media promotes the secretion and deposition of extracellular matrix proteins, such as collagen I, a key protein shown to be upregulated in lung fibrosis.

Figure 2: FMT assay conditions have been optimized to ensure reproducibility and achieve the greatest assay window to allow the study of potential anti-fibrotic therapeutics.


A) 384-well plates were coated with two different coatings and HLF were seeded into wells at an equal density. Cells were fixed and stained using Hoechst 33342 to detect cell nuclei representative images shown of whole well x4 magnification. B) Four assay medias were compared to determine which assay media promotes FMT to provide a robust and significant assay window. HLF were stimulated with TGF- β 1 in four different assay medias. Shown is total integrated intensity normalized to cell number for α -SMA (left) and collagen I (right) of HLF (N=3 donors, each colour point represents a donor) stimulated with 1ng/mL TGF- β 1. Statistical analysis performed; one-way ANOVA with Dunnett's multiple comparisons test compared to Media 1. C) Inclusion of a macromolecular crowding agent (MMC) in the culture media promotes deposition of extracellular matrix and shows an increase in the baseline expression of extracellular collagen I.

TGF- β 1 drives FMT in normal and IPF donors



Inhibition of FMT using ALKi: SB525433*


* Work in development



Newcells Biotech's FMT-assay allows the evaluation of anti-fibrotic therapeutics

Newcells Biotech's FMT assay utilises high-content imaging to determine the effect of therapeutics on disease relevant markers of fibroblast activation (α -SMA) and ECM protein expression (collagen I and EDA-FN1).

Further details available online or contact us at: info@newcellsbiotech.co.uk



NEWCELLS
IN VITRO MODELS TO ACCELERATE CLINICAL TRANSLATION